


Ps™ Trays

Ps™ Instrument Tray



Ps™ Instrument Tray



About Us

PsTwo™

Straightforward and safe

The PsTwo™ System consists of a set of Pedicle Screws, Rods, Set Screws and instruments. PsTwo™ is specially developed for temporary or permanent correction or stabilization of the vertebral column from the thoracic to the sacrum and with the aim of helping consolidation or bone fusion. In combination with a lumbar cage, PsTwo™ bears the intervertebral loads and facilitates osteosynthesis.

PsTwo™ offers the surgeon a straightforward and safe technique. The set offers Mono Screws, Poly Screws and Reduction Screws.

Apollo Swiss Medical

Swiss precision at a reasonable price

Apollo Swiss Medical is a Swiss-based manufacturer of quality spine surgery products. Our experienced medical engineering staff come up with innovative, cost-effective solutions that help surgeons perform implantations with superior effectiveness and reliability.

Our Commitment

Medical practice, therapy center, hospital, academic institution: We care for our customers! To ensure that we can always offer the best possible quality in our products, technology and services, we continually exchange views with practising physicians and researchers and benefit from this wide experience.

Apollo Swiss Medical AG
Spinal Surgery Development & Distribution
Seestrasse 50a
CH-8280 Kreuzlingen

T +41 71 686 89 51
F +41 71 686 89 55

info@apollo-swiss-medical.ch
www.apollo-swiss-medical.ch



Version PS2/01 Stand 7/18



PsTwo™ Spinal Surgery

Pedicle Screw System

Swiss accuracy to satisfy highest requirements



PsTwo™ Implants

The PsTwo™ Implants are produced out of implant grade titanium alloy (Ti6Al4V ELI).
Alle PsTwo™-Schrauben sind kanüliert, fenestriert und haben ein Doppelgewinde für optimierte Fixation.

PsTwo™ Poly Screw



Diameter	Length	Order Number
5,5 mm	35 mm	PT5535-SET
5,5 mm	40 mm	PT5540-SET
5,5 mm	45 mm	PT5545-SET
5,5 mm	50 mm	PT5550-SET
5,5 mm	55 mm	PT5555-SET
6,5 mm	35 mm	PT6535-SET
6,5 mm	40 mm	PT6540-SET
6,5 mm	45 mm	PT6545-SET
6,5 mm	50 mm	PT6550-SET
6,5 mm	55 mm	PT6555-SET
7,5 mm	40 mm	PT7540-SET
7,5 mm	45 mm	PT7545-SET
7,5 mm	50 mm	PT7550-SET
7,5 mm	55 mm	PT7555-SET
8,5 mm	40 mm	PT8540-SINGLE
8,5 mm	45 mm	PT8545-SINGLE
8,5 mm	50 mm	PT8550-SINGLE
8,5 mm	55 mm	PT8555-SINGLE
8,5 mm	60 mm	PT8560-SINGLE
8,5 mm	70 mm	PT8570-SINGLE
8,5 mm	80 mm	PT8580-SINGLE
8,5 mm	90 mm	PT8590-SINGLE

SET beinhaltet jeweils 2 Polyaxialschrauben und 2 Muttern zur Schraube, steril verpackt.
SINGLE Set beinhaltet jeweils 1 Polyaxialschraube und 1 Mutter zur Schraube, steril verpackt.

PsTwo™ Reduction Screw



Diameter	Length	Order Number
5,5 mm	30 mm	RT5530-SET
5,5 mm	35 mm	RT5535-SET
5,5 mm	40 mm	RT5540-SET
5,5 mm	45 mm	RT5545-SET
5,5 mm	50 mm	RT5550-SET
6,5 mm	35 mm	RT6535-SET
6,5 mm	40 mm	RT6540-SET
6,5 mm	45 mm	RT6545-SET
6,5 mm	50 mm	RT6550-SET
6,5 mm	55 mm	RT6555-SET
7,5 mm	40 mm	RT7540-SET
7,5 mm	45 mm	RT7545-SET
7,5 mm	50 mm	RT7550-SET
7,5 mm	55 mm	RT7555-SET
8,5 mm	40 mm	RT8540-SINGLE
8,5 mm	45 mm	RT8545-SINGLE
8,5 mm	50 mm	RT8550-SINGLE
8,5 mm	55 mm	RT8555-SINGLE
8,5 mm	60 mm	RT8560-SINGLE
8,5 mm	70 mm	RT8570-SINGLE
8,5 mm	80 mm	RT8580-SINGLE
8,5 mm	90 mm	RT8590-SINGLE

SET beinhaltet jeweils 2 Repositionsschrauben und 2 Muttern zur Schraube, steril verpackt.
SINGLE Set beinhaltet jeweils 1 Repositionsschraube und 1 Mutter zur Schraube, steril verpackt.

PsOne™ Set Screw



SET beinhaltet jeweils 2 Muttern zur Schraube, steril verpackt.

Order Number

S10-SET


PsOne™ Rod gerade, steriles SET



SET beinhaltet jeweils 2 Stäbe gerade, steril verpackt.
*Stab 480 mm unsteril und einzeln verpackt erhältlich.

Diameter	Length	Order Number
6,0 mm	40 mm	R40-SET
6,0 mm	50 mm	R50-SET
6,0 mm	60 mm	R60-SET
6,0 mm	70 mm	R70-SET
6,0 mm	80 mm	R80-SET
6,0 mm	90 mm	R90-SET
6,0 mm	100 mm	R100-SET
6,0 mm	480 mm	R480*

PsOne™ Rod Standard Curve, steriles SET



SET beinhaltet jeweils 2 Stäbe gebogen, steril verpackt.

Diameter	Length	Order Number
6,0 mm	40 mm	R40SC-SET
6,0 mm	50 mm	R50SC-SET
6,0 mm	60 mm	R60SC-SET
6,0 mm	70 mm	R70SC-SET
6,0 mm	80 mm	R80SC-SET
6,0 mm	90 mm	R90SC-SET
6,0 mm	100 mm	R100SC-SET